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**ULTIMATE
BODY
COMPOSITION
ANALYZER**

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Remove the socks and step on the footplate

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
Wait for the weight measurement finish
the enter your phone number

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
Enter height (cm), years old, gender in order,
then press Enter to go to the next step

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



手臂不要貼著腋下
Your arms must not touch the sides of your body

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
The results will be shown on screen while the test finished and wait for print out

InBody Results	ID	Height	Age	MIF	Exit
Weight		182.5			
Muscle Mass					
Percent Body Fat					
ECW Ratio					
Leg Lean Mass					

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Before using InBody

Few groups of people are not allowed to use InBody :

- Pregnant women
- People with implanted medical electronic devices, such as pacemakers
- People who have joint replacement surgery
- When a child or people with disabilities is using InBody, must be under supervision or assistance.

Before step on InBody, please make sure to remove all heavy items and metal objects on your body, then you can start the measuring.

*Assessment is for reference only not for any diagnostic purpose


InBody 770 Body Composition Analysis Manual

1. Please remove the socks and step on the footplate


2. The screen will show the weight

Wait for the weight measurement finish the enter your phone number


3. Input the personal data

Enter height (cm), years old, gender in order, then press Enter to go to the next step

4. Grip and hold the hand electrodes

Maintain proper posture until the test finish

InBody 770 Body Composition Analysis Manual



5. Perform measurement

Holding the grip electrodes with proper posture


Your arms must not touch the sides of your body

Keep the arms straight


The thighs should not touch

Step on the footplate barefoot

Place the heels on the rear sole electrodes

Place the thumb on the oval electrode and the 4 fingers wrap the surface of the bottom hand electrode

6. The results will be shown on screen while the test finished and wait for print out

InBody Results		ID	Height	Age	MIF	Exit
Weight	73.8 kg		182.5			
Muscle Mass	36.6 kg					
Percent Body Fat	13.1 %					
ECW Ratio	0.368					
Leg Lean Mass	20.2 kg					

Would you like to input ID to save results?

InBody Report Interpretation

Personal information
& Test Date

Body Composition
Analysis

ID	Height	Age	Gender	Test Date / Time	
Jane Doe	163cm	41	Female	2017.03.08. 16 : 47	
Body Composition Analysis					
	Values	Total Body Water	Soft Lean Mass	Fat Free Mass	Weight
1 Total Body Water(L)	35.5 (29.1 ~ 35.5)	35.5	45.6 (37.3 ~ 45.7)	48.3 (39.6 ~ 48.4)	66.4 (48.5 ~ 65.7)
2 Protein (kg)	9.5 (7.8 ~ 9.6)	non-ossious			
3 Minerals (kg)	3.28 (2.69 ~ 3.29)				
4 Body Fat Mass (kg)	18.1 (11.4 ~ 18.3)				

1 Total Body Water

- Body water content

2 Protein

- One of the component build up your muscle

3 Mineral

- Which can maintain body health & metabolism

4 Body Fat Mass

- Body fat weight

5 Soft Lean Mass

- Sum of the body water protein and minerals

6 Fat Free Mass

- Body Weight-Body Fat Mass=Fat Free Mass

7 Weight

- Sum of all values

Values in parentheses are ideal ranges based on your measurements


InBody Report Interpretation

Muscle-Fat Analysis

1. Weight

2. Skeletal Muscle Mass

Muscle responsible for movement and posture ; most easily affected through exercise.


3. Body Fat Mass

The value will divide into three range


Under

Normal

Over

- >> InBody compares your results to other people who have the same [height](#) and [gender](#)
- >> When the value reaches 100%, it reflects your weight is equivalent to a healthy average
- >> If the value reaches 130%, it reflects your weight is **above** the average 30%
- >> In contrast, if the value reaches 70%, it reflects your weight is **below** the average 30%


InBody Report Interpretation

Obesity Analysis

1. BMI

$$\frac{\text{Weight (kg)}}{\text{Height}^2 \text{ (m}^2\text{)}}$$

		Under		Normal			Over					
		10.0	15.0	18.5	21.0	25.0	30.0	35.0	40.0	45.0	50.0	55.0
BMI	Body Mass Index (kg/m ²)	25.0										
PBF	Percent Body Fat (%)	27.2										

- Underweight: <18.5
- Normal: 18.5-22.5
- Risk to overweight: 22.6-24.9
- Overweight: 25-29.9
- Obese: ≥30

2. Percent Body Fat (PBF)

$$\frac{\text{Body fat}}{\text{Mass}} \times 100\%$$

	Values	Total Body Water	Soft Lean Mass	Fat Free Mass	Weight
Total Body Water(L)	35.5 (29.1 - 35.5)	35.5	45.6 (37.3 - 45.7)	48.3 (39.6 - 48.4)	66.4 (48.5 - 65.7)
Protein (kg)	9.5 (7.8 - 9.6)				
Minerals (kg)	3.28 (2.69 - 3.29)				
Body Fat Mass (kg)	18.1 (11.4 - 18.3)				

Body weight

In **Male**, the healthy range of PBF is between 10-20%

In **Female**, the healthy range of PBF is between 18-28%

InBody Report Interpretation

Segmental Lean Analysis

InBody segmental lean analysis provides 5 body parts results :


- ① Right Arm ② Left Arm ③ Trunk ④ Right Leg ⑤ Left Leg

		Based on ideal weight			Based on current weight			ECW Ratio
		Under	Normal	Over	Under	Normal	Over	
Right Arm	(kg)	40	60	80	100	120	140	0.373
	(%)	100			116.6	116.6		
Left Arm	(kg)	40	60	80	100	120	140	0.377
	(%)	107.4			107.4	107.4		
Trunk	(kg)	70	80	90	100	110	120	0.381
	(%)	105.5			105.5	105.5		
Right Leg	(kg)	70	80	90	100	110	120	0.380
	(%)	112.5			112.5	112.5		
Left Leg	(kg)	70	80	90	100	110	120	0.382
	(%)	109.6			109.6	109.6		

- ➔ First values are compared to the people who have the same height. It reflects your muscle weight equivalent to a healthy average when the values equal or above 100%
- ➔ Middle value is the muscle weight of that body part
- ➔ The third values compare your muscle weight to your body weight. It determines if you have enough muscle mass to support your weight. It is enough to support your weight when the values equal or above 100%

InBody Report Interpretation

ECW Ratio Analysis

High SMM/ICW: <0.36

Standard Range: 0.36-0.39

Slight edema: 0.39-0.40

Edema: >0.40

Body Composition History

Graphically tracks the historical trending of your body compositional changes, including :

- ➔ Weight
- ➔ Skeletal Muscle Mass (SMM)
- ➔ Percent Body Fat (PBF)
- ➔ ECW Ratio

InBody Report Interpretation

InBody Score

Total score that reflect the evaluation of body composition. 70-90 score means that the health status is between “good” and “very good” state. A muscular person may score over 100

InBody Score

81 / 100 Points


* Total score that reflects the evaluation of body composition. A muscular person may score over 100 points.

Visceral Fat Area

Visceral fat area determine how much unhealthy visceral fat you have

The ideal range of visceral fat should below 100 VFA cm^2 (Area that highlighted in blue)

Visceral Fat Area

InBody Report Interpretation

Weight Control

Is a recommendation that how to achieve ideal weight after body composition analysis

“+” increase weight

“-” decrease weight

Example : need to lose 3.7kg fat

*InBody do not recommended decrease muscle weight

Weight Control

Target Weight	62.7 kg
Weight Control	- 3.7 kg
Fat Control	- 3.7 kg
Muscle Control	0.0 kg

Segmental Fat Analysis

The values are compared to the people who have the same height

The left values (kg) are the fat mass of this body part

The % values reflect if your fat mass equivalent to a healthy average.

If the value reaches 167%, it reflects your fat mass of this body part is above the average 67%

Segmental Fat Analysis

Right Arm	(1.1kg)	110.9%
Left Arm	(1.2kg)	122.7%
Trunk	(9.0kg)	167.0%
Right Leg	(2.9kg)	119.5%
Left Leg	(2.9kg)	118.0%

Under (較低)

Normal (正常)

Over (較高)

InBody Report Interpretation

Research Parameters

1. Intracellular Water

- Water located within the cell; associated with anabolic processes, overall good cellular health and integrity

2. Extracellular Water

- Water found outside the cell; abnormally high levels often indicative of excess body fat and water retention related to trauma.

3. BMR

- Minimum amount of energy, or calories, required to sustain life at a resting state within a 24-hours period

4. Waist-Hip Rate

- W/H=Waist-Hip Rate
Waist-Hip rate compares your waist measurement to your hip measurement. Higher ratios can mean you have abdominal obesity

5. Body Cell Mass

- Total weight in pounds of all metabolically active tissues of the body

Research Parameters		
1	Intracellular Water	22.0 L (18.0~22.0)
2	Extracellular Water	13.5 L (11.1~13.5)
3	Basal Metabolic Rate	1413 kcal
4	Waist-Hip Ratio	0.83 (0.75~0.85)
5	Body Cell Mass	31.5 kg (25.8~31.6)

InBody Q&A

1. Why the measurement is not successful?

During the measurement, the hand or foot moves away from the electrode contact point. Also, dry skin will cause measure error. It is recommended to thoroughly wipe the palms and the feet with a towelette before taking the measurement.

2. Why does eating affect measurement accuracy?

After a meal, the food and water in the stomach may be considered as fat. The stomach wall interferes with the flow of electricity. Body weight can easily change by 1~2kg due to eating or excreting. Therefore, it is not recommended to test immediately after a meal.

3. Why we cannot exercise before the measurement?

Exercise will change the distribution of water in the body. You can feel that the muscles become stiffness, since the water is concentrated in the muscles. The result may not accurately. Therefore, you should take measurements at least 2~3 hours after exercise.

InBody Q&A

4. Why the result is not the same when test again within 5 minutes ?

Participants may not strictly assume the exact same posture for each test. Moreover, during the rest period, the water distribution of the legs is easy to change, and it is not easy to control the difference. Therefore, it is recommended to wait another 1~2 minutes before taking another test. It is not recommended to test immediately after walking or walking stairs.

5. How often should I take the InBody Test?

When the weight management program is in progress, it is recommended to take another InBody test after 2 weeks. You should have another InBody test within 4 weeks at the latest.

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