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For more than 50 years, spa and salon professionals have used gentle galvanic currents in customized treatment to refresh and energize the skin. These treatments, through a gentle massaging action, help to focus cellular energy and enhance circulation.

Now, small, smart and powerful, the programmable Galvanic Spa™ System II allows you to experience multiple spa benefits at home. Patented technology with self-adjusting galvanic currents and interchangeable conductors, work synergistically with our special formulated products that combine 3-in-1 multi-care to rejuvenate your complexion, revitalize your scalp, and refine your body.

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Revives and invigorates tired
skin within 10 minutes

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
Contains **marigold**, **sea kelp** and **vanilla extracts** that help smooth and condition skin.

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Snap the Face conductor on the Galvanic Spa™ System II.

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Push the selection button once. The "" icon should appear in the display.

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Select a treatment time by pushing the selection button, select the number "2".

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
Contains **arginine** to facilitate recovery from stress while **magnesium** promotes cellular energy.

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Push the selection button twice. The "

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Select a treatment time by pushing the selection button, select the number "3".

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Enjoy 3 minutes spa treatment.

抗皺護理
Anti-Wrinkle
Treatment

撫平幼紋及皺紋
Smooths out the lines and wrinkles

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Features **pro-collagen peptides** that re-activate collagen producing cells, it helps soften moderate to deep lines around mouth, eyes and forehead.

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Clinically proved that the Galvanic Spa™ System II could effectively drive the positively charged vital ingredient - **Peptides** - deeper into the skin, making it more bio-available to the skin to soften lines and wrinkles.

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
Snap the Face
conductor on the
Galvanic Spa™
System II.

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Push the selection
button twice.
The " " icon
should appear in
the display.

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Select a treatment
time by pushing
the selection
button, select the
number "1,2,3 or 5".

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Enjoy 1-5 minutes
spa treatment.

纖體護理
Body Care

擊退橙皮脂，塑造體態美
Repels cellulite, shapes up your body

橙皮脂肪是什麼？ What is cellulite?

- ※ 即「蜂窩組織炎」
- ※ 又稱「橙皮紋」或「茅屋芝士紋」
- ※ 能於身體任何部位出現，尤其皮下脂肪過剩的地方。肌膚表面會不均勻及凹凸不平
- ※ 多出現於大腿上方外側、大腿後側、坐圍、乳房及手臂上方
- ※ 估計約85%女性正受困擾
- ※ Means "cell inflammations"
- ※ Also described as having an "orange peel" or "cottage cheese" appearance
- ※ Can be located anywhere on the body that contains excess subcutaneous fat and appears as uneven, bumpy skin texture
- ※ Commonly seen on the upper outer thighs, the posterior thighs, buttocks, breasts and upper arms
- ※ Estimated that 85% of females are afflicted with the condition

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Contains **theobromine (theobroma cacao extract)**, a molecule shown to target a certain enzyme in fat cells, helping activate fat breakdown; **golden chamomile (chrysanthellum indicum extract)** and testing indicates the **flavonoids** found in golden chamomile help inhibit the storage of fat, increase circulation to remove toxins, even out dimpled areas.

-  分解及減少脂肪囤積
Breakdowns and reduces fat storage
-  加速淋巴管循環，趕走有毒廢物
Increases lymphatic circulation to remove toxins
-  提升肌膚彈性
Increases skin firmness
-  撫平橙皮紋
Evens out dimpled areas

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Clinical study shows that capillary and lymphatic circulation improvement can also assist with the treatment of cellulite. An extract of **chrysanthellum indicum** can reduce edema and inflammation and promote tissue drainage and the elimination of toxins.

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Snap the Body conductor on the Galvanic Spa™ System II.

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大腿尺寸減少 Reduces in thigh diameter	73%
橙皮脂肪明顯改善 Visible improvement in appearance of cellulite	50%
整體橙皮紋贅肉改善情況 Improvement in overall cellulite condition	46%
肌膚平滑度 Exhibits a smoother skin surface	35%

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Used in conjunction with Dermatic Effects™ Body Contouring Lotion, it almost **doubles** the reduction in the extent of the cellulite condition.

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Contains **hibiscus** and **malvaceae extracts** that help inhibit fat production and stimulate fat breakdown. It also helps prevent the appearance of cellulite and skin aging.

專業秀髮護理
Professional Hair
Treatment


讓您擁有健康豐盈秀髮，
散發自信光采
Helps your hair look healthy and abundant,
shows off your confidence in no doubt

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Each hair follicle cycles through 3 phases:

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Clinical study demonstrated that if we could provide abundant nutrient and sustenance to the hair bulb, it helps extent the growth phase in the hair life cycle.

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Formulated with patented technology **Tricalgoxyl®**, not only remineralizes the scalp, but also increases the availability of key nutrients and vitamins to the hair follicle. It helps revitalizes and strengthens hair in improving volume and shine as well.

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Snap the Scalp conductor on the Galvanic Spa™ System II.

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A U.S. Corporation with Limited Liability

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All Right Reserved

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Printed in Hong Kong

04133150

May 2008