

# ageLOC<sup>®</sup> boost<sup>™</sup> SYSTEM



Let your skin radiate brilliance with an innovative device system from Nu Skin. ageLOC Boost system features unique, variable pulse microcurrent technology and innovative ingredients to visibly brighten and awaken your skin. Now that's simply brilliant!

## MEET AGELOC BOOST SYSTEM

Say goodbye to dull, tired-looking skin and hello to skin that is bright, bouncy, brilliant. Introducing ageLOC Boost, a breakthrough beauty device that works with the ageLOC Boost Activating Treatment to bring out your skin's natural beauty. We're talking visibly brighter, plumper, bouncier skin with youthful radiance and a more even-toned complexion. Effortless to use (just two minutes a day), ageLOC Boost system gives you a beauty boost with every use.

## WHY YOU'LL LOVE IT

- Boosts your skin's natural beauty for brighter, plumper, and bouncier skin.
- Potentially reduces the appearance of dark spots and hyperpigmentation.
- May improve and even skin tone for a brighter, more radiant complexion.
- May help in the maintenance of youthful-looking skin, easing the signs of premature aging.
- Hydrates, so skin feels softer and appears fuller.

## WHY YOU'LL LOVE IT EVEN MORE

- Features patent-pending Variable Pulse Microcurrent Technology.
- Delivers lasting benefits.<sup>1</sup>
- Keeps you at your selfie- and video-chat-ready best!
- Super easy and quick to use (just two minutes a day).
- Effortless add-on to most skin care routines.
- Features a bottle that is proudly made with 15% post-consumer resin (PCR).

## WHAT POWERS IT

- **Spotless Concentrate**—may reduce the appearance of hyperpigmentation and dark spots for a more even-looking skin tone. It contains Sea Daffodil Extract, Vitamin C, Vitamin B3 and Niacinamide.
- **Plumping Blend**—hydrates and softens, while giving skin a fuller, plumper appearance. It contains Watermelon Rind, Apple Starch, Lentil Extract and Hyaluronic Acid in 3 different weights of Sodium Hyaluronate and another form as Sodium Acetylated Hyaluronate.
- **ageLOC<sup>®</sup> Blend**—Nu Skin's exclusive technology that targets signs of aging to help maintain a youthful, vibrant appearance.

## WHAT YOU'LL RECEIVE

- ageLOC Boost device with magnetic charger
- ageLOC Boost Activating Treatment

## WHAT'S VARIABLE PULSE MICROCURRENT TECHNOLOGY?

The ageLOC Boost device sends intermittent, variable pulses to the skin, eight times a second. The pulses are charge-balanced, so they switch from positive to negative every 2.4 seconds. When used with ageLOC Boost Activating Treatment, the microcurrent optimizes and activates this formula to its maximum potential.

## AGELOC BOOST ACTIVATING TREATMENT—ADDING OOMPH TO YOUR SKIN CARE

Formulated to work with the device, ageLOC Boost Activating Treatment will help you achieve a glowing, youthful complexion. It contains efficacious brightening

<sup>1</sup>With regular use of product

# ageLOC® boost™

ingredients that target the appearance of hyperpigmentation, a hydrating plumping blend, and our proprietary ageLOC Blend. The outcome? Visibly bouncier skin with a healthy radiance.

## DEMONSTRATED RESULTS

Nu Skin Research and Development commissioned a third-party clinical research laboratory to identify the visible skin benefits of using ageLOC Boost with ageLOC Boost Activating Treatment. The 12-week study included Asian females between the ages of 25 and 40 years who used it once a day.

### Clinical Evaluation

A clinical grader noted statistically significant results of:

85% improvement in radiance

85% improvement in skin tone evenness

78% improvement in firmness

77% improvement in bounce

75% improvement in texture

73% improvement in plumpness

### Self-Perception

Subjects noted statistically significant changes in their skin's appearance by week 1 and improvements continued to grow throughout the study.

At the end of the study, they self-evaluated with the following results:

108% improvement in appearance of dark spots

106% improvement in hydration

103% improvement in pore visibility

97% improvement in smoothness

95% improvement in radiance

93% improvement in bounce

93% improvement in plumpness

100% of subjects also agree with the following statements:

"My skin appears brighter"

"My skin looks dewy"

"My skin looks baby plump"

"My skin looks 'lit from within'"


## HOW TO USE IT

Super easy to use, ageLOC Boost system transforms your skin from blah to brilliant. After cleansing and toning, follow these easy steps:

1. Turn on device by swiping finger up full length of touch surface.



2. Dispense three pumps of ageLOC Boost Activating Treatment onto the device treatment surface. Spread over half of the face.



3. Begin the 1 minute treatment by swiping finger up full length of touch surface. Glide treatment surface over the half of your face in an upward/outward motion.



4. Repeat STEP 2 & STEP 3 on the other half of the face for another 1 minute.



5. Massage excess ageLOC Boost Activating Treatment into skin. (Device will automatically turn off)

And that's it—you're done! Follow up with your favorite serum, moisturizer, and sunscreen. You'll immediately notice your skin feels more hydrated and softer.

Please note: only use ageLOC Boost Activating Treatment with the ageLOC Boost device, since it was scientifically formulated to conduct the unique variable pulse microcurrent.

## DID YOU KNOW

ageLOC Boost is an award-winning product recognized on the global landscape by prestigious organizations:

- Gold Winner for Golden Bridge Awards 2020 in Consumer Products
- Silver Winner for Consumer World Awards 2020 in Best Consumer Product
- Bronze 2020 Stevie Winner for Best New Product or Service of the Year in Consumer Products, Household Products
- Winner of BIG Awards for Business 2020 in New Product of the Year



### TRY IT WITH THESE

- **ageLOC LumiSpa®**—a first-of-its-kind, dual-action skin care system that helps to deliver dramatic skin renewal and deep cleansing in 1 Step for 2 Minutes with 7 Benefits.
- **ageLOC Galvanic Spa®**—a single device designed for multiple uses depending on the compatible products you pair with. It uses direct/galvanic currents to help visibly rejuvenate key areas of your face, body and scalp.
- **ageLOC Me®**—creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five uniquely formulated anti-aging products, custom selected from almost 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.

### SMART DESIGN

To provide a positive user experience, we created ageLOC Boost in partnership with an award-winning, world-renowned design firm. The device features an ergonomic feel that's comfortable and easy to use.

### FREQUENTLY ASKED QUESTIONS

#### When should I use ageLOC Boost?

You can use ageLOC Boost anytime. However, we recommend using it as part of your morning skin care routine, so you'll get a beauty boost that lasts throughout the day.

#### How are ageLOC Boost and ageLOC Galvanic Spa® different?

You can use ageLOC Boost every day, whereas you should use ageLOC Galvanic Spa two to three times a week. The main difference between the two devices is they use a different type of current. ageLOC Boost uses exclusive variable pulse microcurrents to optimize skin, which helps reveal a more youthful complexion. ageLOC Galvanic Spa uses direct microcurrents to deliver treatment products to visibly improve aged skin.

#### Can I use ageLOC Boost in the shower?

You should not use ageLOC Boost in the shower or bath because the device is not waterproof (it's water-resistant). Also, you'll want to use the system after showering because Activating Treatment is a leave-on product, not one you wash off.

#### How does ageLOC Boost support Nu Skin's commitment to become more sustainable?

ageLOC Boost's Activating Treatment bottle is made with 15% post-consumer resin (PCR), and the product and device cartons are made from virgin paper. Please check your local capabilities and recycle the bottle and cartons where possible. By purchasing ageLOC Boost, you're joining us in taking an active, engaged stride to ensure the future is bright for both the planet and each of us. Please visit [www.nuskin.com/sustainability](http://www.nuskin.com/sustainability) for more information and details on how to recycle.

#### Can I use any product with ageLOC Boost?

To enjoy the best treatment experience, you should only use ageLOC Boost Activating Treatment with your device. That's because we scientifically formulated Activating Treatment to optimally conduct the unique variable pulse microcurrent.

#### Why do I need to put ageLOC Boost Activating Treatment on the device's treatment head, rather than directly on my skin?

Dry hands are a must when using the device, which is why we've developed a treatment protocol that lets you dispense ageLOC Boost Activating Treatment onto the device head, so your hands remain dry.

### INGREDIENTS

ageLOC Boost Activating Treatment  
Water (Aqua), Glycerin, Butylene Glycol, Dimethicone, Niacinamide, Tetrahexyldecyl Ascorbate, Polyglyceryl-6 Distearate, Sodium Acetylated Hyaluronate, Sodium Hyaluronate, Sodium PCA, Citrullus Lanatus (Watermelon) Fruit Extract, Lens Esculenta (Lentil) Fruit Extract, Pyrus Malus (Apple) Fruit Extract, Pancreaticum Maritimum Extract, Sodium Lactate, Ergothioneine, Tocopheryl Acetate, Jojoba Esters, Cetyl Alcohol, Lecithin, Mannitol, Polyglyceryl-6 Polyricinoleate, Polyglyceryl-3 Beeswax, Xanthan Gum, Hydroxypropyl Methylcellulose Stearoxy Ether, Sodium Acrylates Copolymer, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Propanediol, Chlorphenesin, Ethylhexylglycerin, Phenoxyethanol, Benzoic Acid, Tetrasodium Glutamate Diacetate, Aminomethyl Propanol.

